



MICROPLATES

AFI



AFI-RBHDWN-150 : 8 MICROPLATES

Model	V	R
Max Capacity	-	-
Angle	90°	
RPM Max	4 500	4 500
ACR Max	2 400	2 400



Microplates

Capacity (ml)	-	-	-	-
Max Height (mm)	-	-	55	55
Tubes per adapter / rotor	-	-	4 8	1 2

Accessories



Hole \varnothing x L (mm)	-	-	-
Radius	106	-	106
Max Charge (gr)	-	-	575
Ref.	AFI-RBHDW-N150	AFI-RHDW-N150	AFI-BHDW-N150*

* Set of 2