



# ACHIEVE HIGH SURVIVAL RATES WITH The Cryotec Method

Cryotec products are designed for safe and effective operability



## Solutions enable safe recovery

- Safe and long-lasting due to solutions being free from HSA and proteins
- Low impact on cells due to the use of endotoxin-free trehalose and a low concentration of CPA
- The specific gravity of each solution makes it easy to identify the completion of the equilibration in vitrification steps
- Hydroxypropyl cellulose (HPC) increases viscosity and slows down the movement of embryos



The Cryotec Method is a **safe and effective** technique for vitrification and warming of oocytes and embryos. It offers a unified protocol for oocytes and embryos of all stages leading to **high survival rates**.



**NO TIME PRESSURE**



**ONE PROTOCOL FOR ALL**



**GENTLE ON OOCYTES AND EMBRYOS**



**IMPROVED OPERABILITY**

## Our Products

### Vitrification Kit 101 (For 3 applications)



1.0 ml Equilibration Solution (ES) x1  
1.0 ml Vitrification Solution (VS) x2  
Cryotec x4  
Vitri Plates x3

### Vitrification Solution Set 110 (For 10 applications)



1.8 ml Equilibration Solution (ES) x2  
1.8 ml Vitrification Solution (VS) x4

### Vitri Plate (Set of 10)



### Warming Kit 102 (For 1 application)



1.8 ml Thawing Solution (TS) x1  
0.5 ml Diluent Solution (DS) x1  
1.0 ml Washing Solution (WS) x1  
Warm Plate x1

### Warming Solutions Set 205 (For 5 applications)



1.8 ml Thawing Solution (TS) x5  
1.8 ml Diluent Solution (DS) x1  
1.8 ml Washing Solution (WS) x2

### Warm Plate (Set of 10)



### Cryotec

10 Cryotec vitrification devices



Available in blue, green, yellow, white and red

## How to order

Please contact us by email [globalsales@reprolife.jp](mailto:globalsales@reprolife.jp) or scan QR code to visit our contact page



**Note:** Cryotec products are designed for simple, safe and effective operability to give the best chance for successful vitrification and warming of oocytes and embryos. REPROLIFE Inc. is ISO 13485 certified. We only ship products that have passed a variety of rigorous tests.

This brochure was updated in July 2023

Produced by REPROLIFE Inc.



Disclaimer: Some products are not available in the USA, EU and other countries. Contact us to inquire.

